




Floris

Floris IS A MANUFACTURED FILLER USING HYALURONICACID (HA) DERIVED FROM FREE ANIMAL SUBSTANCES WITH HIGH PURITYTOHAVEBIOCOMPATIBILITY IN HUMAN TISSUE.

HYALURONIC ACID IS A COMPONENT THAT HELPS SUPPLY WATER AND REGENERATE CELLS IN THE SKIN. FLORIS IS COMPOSED OF RETI CULAR DERIVATIVES OF SWELLING DEGREES THAT CONSTANTLY MAINTAIN ELASTICITY INSIDE THE BODY.

"DIFFERENTIATED MANUFACTURING TECHNOLOGY" PREMIUM HYALURONIC ACID, INTERNATIONALLY RECOGNIZED BY FDA AND EDQM, IS BEING USED.



Floris Core technology

Definition of MSFT (Multiple Swelling Formulation Technology)

01


MSFT is a method of blending multiple swelling ratio, and it maintains excellent volume effects while ensuring safety by minimizing crosslinking agent.



Minimal use of crosslinking agent





Increase hydrophilicity by blending gels with different swelling ratio



Safely maintain volume by minimizing crosslinking agent and increasing hydrophilicity

02

Floris, produced by MSFT technology, has tissue adhesion and excellent volume effect.



03

Floris is an innovative product that minimizes toxicity and increases safety with the use of its unique technology MSFT


	Type	FLORIS	R filler	J filler
Toxicity	Endotoxin_Pyrogen	Under 0.05EU/mL	Under 0.5EU/mL	Under 0.5EU/mL
	Protein Load_Foreign body reaction	Under 8µg/mL	Under 17~13µg/mL	Under 39µg/mL
	BDDE_residue	Under 2ppm	Under 2ppm	Under 2ppm

Floris minimizes free radical formation by significantly reducing proteins that may cause hypersensitivity reactions



→ *Floris*'s protein content: under 8 µg /mL

Floris has a low unit level of endotoxin which holds an outstanding safety profile.


→ *Floris*'s endotoxin unit level: under 0.05EU/mL



04 *Floris*, produced by MSFT technology, has tissue adhesion and excellent volume effect.



05 *Floris* results in a superior outcome in Syringeability compared to other companies.



06 *Floris* combines all the necessary qualities: Safety, Efficacy, Sustainability.

Multiple Swelling Formulation Technology

**Soft
injection
force**

**Excellent
volume
effect**

**In-vivo
safety**

**Superior
adhesion**

**Outstanding
sustainability**

**Toxicity
minimization**

**Hyaluronic
acid
derivative**



Specification

	Floris fine	Floris deep	Floris ultra	Floris S	Floris M	Floris H
Composition	HA 20mg/mL	HA 20mg/mL	HA 20mg/mL	HA 20mg/mL Lidocaine %0.3	HA 24mg/mL Lidocaine %0.3	HA 24mg/mL Lidocaine %0.3
Indications	Superficial lines : Worry lines, Crow's feet, Neck Winkles	Deep nasolabial folds, Face contouring and augmentation	Severely deep lines, Face contouring and augmentation	Superficial lines : Worry lines, Crow's feet, Neck Winkles	Deep nasolabial folds, Face contouring and augmentation	Severely deep lines, Face contouring and augmentation
Injection Depth	Superficial dermis	Mid/Deep dermis, Subcutaneous	Deep dermis, Subcutaneous	Superficial dermis	Mid/Deep dermis, Subcutaneous	Deep dermis, Subcutaneous
Volume	1.1mL	1.1mL	1.1mL	10mL	10mL	10mL
Needle Size	30 G	27 G	27 G	30 G	27 G	27 G
Duration	12~8 month	12~8 month	18~10 month	12~8 month	12~8 month	18~10 month

